

EDF-4096

Modified Method 8280 Matrix Spiking Solution
(All 17 toxic congeners)1.2 mL
in Nonane

Unlabeled Dioxins & Furans	Concentration (ng/μl)
2,3,7,8-TCDD	2.5
2,3,7,8-TCDF	2.5
1,2,3,7,8-PeCDD	6.25
1,2,3,7,8-PeCDF	6.25
2,3,4,7,8-PeCDF	6.25
1,2,3,4,7,8-HxCDD	6.25
1,2,3,4,7,8-HxCDF	6.25
1,2,3,6,7,8-HxCDD	6.25
1,2,3,6,7,8-HxCDF	6.25
1,2,3,7,8,9-HxCDD	6.25
1,2,3,7,8,9-HxCDF	6.25
2,3,4,6,7,8-HxCDF	6.25
1,2,3,4,6,7,8-HpCDD	6.25
1,2,3,4,6,7,8-HpCDF	6.25
1,2,3,4,7,8,9-HpCDF	6.25
OCDD	12.5
OCDF	12.5