



Cambridge Isotope Laboratories, Inc.

3 Highwood Drive, Tewksbury, MA 01876 USA

PH: 800.322.1174 (N. America) PH: 978-749-8000

Fax: 978-749-2768 WEB: www.isotope.com

EC-4133

DSJ PCB Mixture

1.2 mL in isoctane

Unlabeled	IUPAC	Concentration (ng/mL)
2,2',4-TriCB	17	500
2,2',5-TriCB	18	2000
2,4,4'-TriCB	28	2000
2,4',5-TriCB	31	1500
2',3,4-TriCB	33	2000
2,2',3,5'-TetraCB	44	2000
2,2',4,5'-TetraCB	49	2000
2,2',5,5'-TetraCB	52	2000
2,3',4',5-TetraCB	70	2000
2,4,4',5-TetraCB	74	2000
2,2',3,3',4-PentaCB	82	500
2,2',3,4,5'-PentaCB	87	2000
2,2',3,5',6-PentaCB	95	1000
2,2',4,4',5-PentaCB	99	2000
2,2',4,5,5'-PentaCB	101	2000
2,3,3',4,4'-PentaCB	105	500
2,3,3',4',6-PentaCB	110	2000
2,3',4,4',5-PentaCB	118	2000
2,2',3,3',4,4'-HexaCB	128	2000
2,2',3,3',4,6'-HexaCB	132	1000
2,2',3,4,4',5'-HexaCB	138	2000
2,2',3,4',5',6-HexaCB	149	2000
2,2',3,5,5',6-HexaCB	151	2000
2,2',4,4',5,5'-HexaCB	153	2000
2,3,3',4,4',5-HexaCB	156	2000
2,3,3',4,4',6-HexaCB	158	500
3,3',4,4',5,5'-HexaCB	169	2000
2,2',3,3',4,4',5-HeptaCB	170	2000
2,2',3,3',4,4',6-HeptaCB	171	2000

2,2',3,3',4',5,6-HeptaCB	177	2000
2,2',3,4,4',5,5'-HeptaCB	180	2000
2,2',3,4,4',5',6-HeptaCB	183	2000
2,2',3,4',5,5',6-HeptaCB	187	2000
2,3,3',4,4',5',6-HeptaCB	191	2000
2,2',3,3',4,4',5,5'-OctaCB	194	2000
2,2',3,3',4,4',5,6-OctaCB	195	2000
2,2',3,3',4,5,5',6'-OctaCB	199	1500
2,3,3',4,4',5,5',6-OctaCB	205	2000
2,2',3,3',4,4',5,5',6-NonaCB	206	2000
2,2',3,3',4,5,5',6,6'-NonaCB	208	2000
DecaCB	209	2000