

EDF-5008-50

US EPA Method 8290 Matrix Spiking Solution
unlabeled (1:50 dilution of EDF-5008)0.1 mL
in Nonane

Labeled Component	Concentration (ng/mL)
2,3,7,8-Tetrachlorodibenzo-p-dioxin	2
2,3,7,8-Tetrachlorodibenzofuran	2
1,2,3,7,8-Pentachlorodibenzo-p-dioxin	5
1,2,3,7,8-Pentachlorodibenzofuran	5
2,3,4,7,8-Pentachlorodibenzofuran	5
1,2,3,4,7,8-Hexachlorodibenzo-p-dioxin	5
1,2,3,4,7,8-Hexachlorodibenzofuran	5
1,2,3,6,7,8-Hexachlorodibenzo-p-dioxin	5
1,2,3,6,7,8-Hexachlorodibenzofuran	5
1,2,3,7,8,9-Hexachlorodibenzo-p-dioxin	5
1,2,3,7,8,9-Hexachlorodibenzofuran	5
2,3,4,6,7,8-Hexachlorodibenzofuran	5
1,2,3,4,6,7,8-Heptachlorodibenzo-p-dioxin	5
1,2,3,4,6,7,8-Heptachlorodibenzofuran	5
1,2,3,4,7,8,9-Heptachlorodibenzofuran	5
Octachlorodibenzo-p-dioxin	10
Octachlorodibenzofuran	10