



Cambridge Isotope Laboratories, Inc.

3 Highwood Drive, Tewksbury, MA 01876 USA

PH: 800.322.1174 (N. America) PH: 978-749-8000

Fax: 978-749-2768 WEB: www.isotope.com

EDF-4139-10

Extraction Standard for EN-1948-2

2 x 5 mL in nonane

Labeled	Concentration (ng/mL)
2,3,7,8-TetraCDD ($^{13}\text{C}_{12}$, 99%)	100
2,3,7,8-TetraCDF ($^{13}\text{C}_{12}$, 99%)	100
1,2,3,7,8-PentaCDD ($^{13}\text{C}_{12}$, 99%)	100
2,3,4,7,8-PentaCDF ($^{13}\text{C}_{12}$, 99%)	100
1,2,3,4,7,8-HexaCDD ($^{13}\text{C}_{12}$, 99%)	100
1,2,3,4,7,8-HexaCDF ($^{13}\text{C}_{12}$, 99%)	100
1,2,3,6,7,8-HexaCDD ($^{13}\text{C}_{12}$, 99%)	100
1,2,3,6,7,8-HexaCDF ($^{13}\text{C}_{12}$, 99%)	100
2,3,4,6,7,8-HexaCDF ($^{13}\text{C}_{12}$, 99%)	100
1,2,3,4,6,7,8-HeptaCDD ($^{13}\text{C}_{12}$, 99%)	200
1,2,3,4,6,7,8-HeptaCDF ($^{13}\text{C}_{12}$, 99%)	200
OctaCDD ($^{13}\text{C}_{12}$, 99%)	200
OctaCDF ($^{13}\text{C}_{12}$, 99%)	200