



Cambridge Isotope Laboratories, Inc.

3 Highwood Drive, Tewksbury, MA 01876 USA

PH: 800.322.1174 (N. America) PH: 978-749-8000

Fax: 978-749-2768 WEB: www.isotope.com

ES-5611

Method 537 PFC Surrogate Standard Mix ( $^{13}\text{C}$ ; D)

1.2 mL in methanol

<b>Labeled</b>	<b>Concentration (ng/mL)*</b>
Perfluorohexanoic acid (PFHxA), sodium salt ( $^{13}\text{C}_6$ , 99%)	1000
Perfluorodecanoic acid (PFDA) ( $^{13}\text{C}_9$ , 99%)	1000
<i>N</i> -Ethylperfluorooctanesulfonamidoacetic acid ( <i>N</i> -EtFOSAA) ( <i>N</i> -ethyl- $\text{D}_5$ , 98%)	4000

**\*Gravimetric Concentration is based on the salt form, where applicable**