



Cambridge Isotope Laboratories, Inc.

3 Highwood Drive, Tewksbury, MA 01876 USA

PH: 800.322.1174 (N. America) PH: 978-749-8000

Fax: 978-749-2768 WEB: www.isotope.com

EC-5085

Toxic and Predominant PCB PAR Solution

1.2 mL in nonane

Unlabeled	IUPAC	Concentration (ng/mL)
2,2',5-TriCB	18	250
2,4,4'-TriCB	28	250
2,2',3,5'-TetraCB	44	250
2,2',4,5'-TetraCB	49	250
2,2',5,5'-TetraCB	52	250
2,3',4,4'-TetraCB	66	250
2,4,4',5-TetraCB	74	250
2,2',3,4,5'-PentaCB	87	250
2,2',4,4',5-PentaCB	99	250
2,2',4,5,5'-PentaCB	101	250
2,3,3',4,4'-PentaCB	105	250
2,3,3',4',6-PentaCB	110	250
2,3',4,4',5-PentaCB	118	250
2,2',3,3',4,4'-HexaCB	128	250
2,2',3,4,4',5-HexaCB	138	250
2,2',3,4',5,5'-HexaCB	146	250
2,2',3,4',5',6-HexaCB	149	250
2,2',3,5,5',6-HexaCB	151	250
2,2',4,4',5,5'-HexaCB	153	250
2,3,3',4,4',5-HexaCB	156	250
2,3,3',4,4',5'-HexaCB	157	250
2,3,3',4,4',6-HexaCB	158	250
2,3',4,4',5,5'-HexaCB	167	250
2,2',3,3',4,4',5-HeptaCB	170	250
2,2',3,3',4,5,5'-HeptaCB	172	250
2,2',3,3',4',5,6-HeptaCB	177	250
2,2',3,3',5,5',6-HeptaCB	178	250
2,2',3,4,4',5,5'-HeptaCB	180	250
2,2',3,4,4',5',6-HeptaCB	183	250

2,2',3,4',5,5',6-HeptaCB	187	250
2,3,3',4,4',5,5'-HeptaCB	189	250
2,2',3,3',4,4',5,5'-OctaCB	194	250
2,2',3,3',4,4',5,6-OctaCB	195	250
2,2',3,3',4,4',5',6-OctaCB	196	250
2,2',3,3',4,5,5',6'-OctaCB	199	250
2,2',3,4,4',5,5',6-OctaCB	203	250
2,2',3,3',4,4',5,5',6-NonaCB	206	250
DecaCB	209	250