

ES-5529

Phosphorus Flame Retardant Solution (D, 98%)

1.2 mL  
in Acetonitrile

Labeled Compound	Concentration (µg/mL)
Tris(2-chloroethyl) phosphate (D <sub>12</sub> , 98%)	10
Tris(2-chloroisopropyl) phosphate (D <sub>18</sub> , 98%)	10
Tris(1,3-dichloro-2-propyl) phosphate (D <sub>15</sub> , 98%)	10
Triethyl phosphate (D <sub>15</sub> , 98%)	10
Tripropyl phosphate (D <sub>21</sub> , 98%)	10
Tributyl phosphate (D <sub>27</sub> , 98%)	10
Triphenyl phosphate (D <sub>15</sub> , 98%)	10