

EC-4979

Method 1668A Labeled Injection Internal Standard Solution

1.2 mL
in Nonane

Labeled PCBs	IUPAC*	Concentration ($\mu\text{g/mL}$)
2,5-DiCB ($^{13}\text{C}_{12}$,99%)	9L	5.0
2,2',5,5'-TetraCB ($^{13}\text{C}_{12}$,99%)	52L	5.0
2,2',4,5,5'-PentaCB ($^{13}\text{C}_{12}$,99%)	101L	5.0
2,2',3,4,4',5'-HexaCB ($^{13}\text{C}_{12}$,99%)	138L	5.0
2,2',3,3',4,4',5,5'-OctaCB ($^{13}\text{C}_{12}$,99%)	194L	5.0