

EC-5397

Co-PCB Syringe Standard ($^{13}\text{C}_{12}$, 99%)1.2 mL
in Nonane

Labeled Compounds	IUPAC#	Concentration (ng/mL)
$^{13}\text{C}_{12}$ -2,3',4',5-Tetrachlorobiphenyl	70	20
$^{13}\text{C}_{12}$ -2,3,3',5,5'-Pentachlorobiphenyl	111	20
$^{13}\text{C}_{12}$ -2,2',3,4,4',5'-Hexachlorobiphenyl	138	20
$^{13}\text{C}_{12}$ -2,2',3,3',5,5',6-Heptachlorobiphenyl	178	20