

ECN-5489-CS4 PCN Calibration Solution CS4
(unlabeled/¹³C₁₀, 99%)

0.5 mL
in Nonane

Unlabeled Component	Concentration (ng/mL)
1,2,3,4-TetraCN	2
1,2,3,5,7-PentaCN	2
1,2,3,4,6,7-HexaCN	2
1,2,3,5,6,7-HexaCN	2
1,2,3,5,6,8-HexaCN	2
1,2,3,4,5,6,7-HeptaCN	2
OctaCN	2
Labeled Component	
1,2,3,4-TetraCN (¹³ C ₁₀ , 99%)	10
1,2,3,5,7-PentaCN (¹³ C ₁₀ , 99%)	10
1,2,3,4,5,7-HexaCN (¹³ C ₁₀ , 99%)	10
1,2,3,5,6,7-HexaCN (¹³ C ₁₀ , 99%)	10
1,2,3,4,5,6,7-HeptaCN (¹³ C ₁₀ , 99%)	10
OctaCN (¹³ C ₁₀ , 99%)	10
1,2,3,4-TCDD (¹³ C ₆ , 99%)	25