

ECN-5497

PCN Native PAR Solution

1.2 mL  
in Nonane

Unlabeled Component	PCN	Concentration (ng/ml)
1,2,3,4-TetraCN	27	1000
1,2,3,5,7-PentaCN	52	1000
1,2,3,4,6,7-HexaCN	66	1000
1,2,3,5,6,7-HexaCN	67	1000
1,2,3,5,6,8-HexaCN	68	1000
1,2,3,4,5,6,7-HeptaCN	73	1000
OctaCN	75	1000