

ES-5484

CDC OH-PAH Native PAR Standard (unlabeled)

1.2 mL
in Toluene

Unlabeled Component	Concentration (ng/mL)
1-Hydroxynaphthalene	4000
2-Hydroxynaphthalene	4000
2-Hydroxyfluorene	1000
3-Hydroxyfluorene	1000
9-Hydroxyfluorene	1000
1-Hydroxyphenanthrene	1000
2-Hydroxyphenanthrene	1000
3-Hydroxyphenanthrene	1000
4-Hydroxyphenanthrene	1000
1-Hydroxypyrene	1000