

ES-5543

US EPA Method 8276 Toxaphene Composite Stock Standard

1.2 mL
in Nonane

| Unlabeled Component | Concentration (ng/mL) |
|---------------------|-----------------------|
| Hx-Sed | 1000 |
| Hp-Sed | 1000 |
| Parlar 26 | 1000 |
| Parlar 40 | 1000 |
| Parlar 41 | 1000 |
| Parlar 44 | 1000 |
| Parlar 50 | 1000 |
| Parlar 62 | 1000 |