

EC-5000

WHO Coplanar &amp; Mono-Ortho PCBs &amp; 170/180

1.2 mL  
in Isooctane

<b>Unlabeled PCBs</b>	<b>IUPAC*</b>	<b>Concentration (<math>\mu\text{g/mL}</math>)</b>
3,3',4,4'-TetraCB	77	2.0
3,4,4',5-TetraCB	81	2.0
2,3,3',4,4'-PentaCB	105	2.0
2,3,4,4',5-PentaCB	114	2.0
2,3',4,4',5-PentaCB	118	2.0
2',3,4,4',5-PentaCB	123	2.0
3,3',4,4',5-PentaCB	126	2.0
2,3,3',4,4',5-HexaCB	156	2.0
2,3,3',4,4',5'-HexaCB	157	2.0
2,3',4,4',5,5'-HexaCB	167	2.0
3,3',4,4',5,5'-HexaCB	169	2.0
2,2',3,3',4,4',5-HeptaCB	170	2.0
2,2',3,4,4',5,5'-HeptaCB	180	2.0
2,3,3',4,4',5,5'-HeptaCB	189	2.0