

EC-5371

EN-1948-4 PCB Recovery Standard ($^{13}\text{C}_{12}$, 99%)1.2 mL
in Nonane

| ^{13}C -Labeled Compounds | IUPAC# | Concentration (ng/mL) |
|--|--------|-----------------------|
| $^{13}\text{C}_{12}$ -2,3',4',5-TetraCB | 70 | 100 |
| $^{13}\text{C}_{12}$ -2,3,3',5,5'-PentaCB | 111 | 100 |
| $^{13}\text{C}_{12}$ -2,2',3,3',4,4',5-HeptaCB | 170 | 100 |