

EC-5495

Marker-7 PCB Mixture (with PCB-118) unlabeled

1.2 mL
in Isooctane

Unlabeled Compounds	IUPAC#	Concentration (ng/mL)
2,4,4'-TriCB	28	1000
2,2',5,5'-TetraCB	52	1000
2,2',4,5,5'-PentaCB	101	1000
2,3',4,4',5-PentaCB	118	1000
2,2',4,4',5,5'-HexaCB	153	1000
2,2',3,4,4',5'-HexaCB	138	1000
2,2',3,4,4',5,5'-HeptaCB	180	1000