

EC-5601

EPA Method 1668A/B/C Alt. Injection
Internal Standard Solution1.2 mL
in Isooctane

Unlabeled	IUPAC	Concentration (ng/mL)
2,5-DiCB	9	1000
2,4',5-TriCB	31	1000
2,4',6-TriCB	32	1000
2,2',5,5'-TetraCB	52	1000
2,2',4,5,5'-PentaCB	101	1000
3,3',4,5,5'-PentaCB	127	1000
2,2',3,4,4',5'-HexaCB	138	1000
2,2',3,4,4',5,5'-HeptaCB	180	1000
2,2',3,3',4,4',5,5'-OctaCB	194	1000