

ECN-5489-CS2 PCN Calibration Solution CS2
(unlabeled/¹³C₁₀, 99%)

0.5 mL
in Nonane

Unlabeled Component	Concentration (ng/mL)
1,2,3,4-TetraCN	0.2
1,2,3,5,7-PentaCN	0.2
1,2,3,4,6,7-HexaCN	0.2
1,2,3,5,6,7-HexaCN	0.2
1,2,3,5,6,8-HexaCN	0.2
1,2,3,4,5,6,7-HeptaCN	0.2
OctaCN	0.2
Labeled Component	
1,2,3,4-TetraCN (¹³ C ₁₀ , 99%)	10
1,2,3,5,7-PentaCN (¹³ C ₁₀ , 99%)	10
1,2,3,4,5,7-HexaCN (¹³ C ₁₀ , 99%)	10
1,2,3,5,6,7-HexaCN (¹³ C ₁₀ , 99%)	10
1,2,3,4,5,6,7-HeptaCN (¹³ C ₁₀ , 99%)	10
OctaCN (¹³ C ₁₀ , 99%)	10
1,2,3,4-TCDD (¹³ C ₆ , 99%)	25