

EM-4173

Method 1653A Labeled Chlorophenolic Derivatives Mixture

2 x 1 mL

(Note: Unlabeled Internal Standard, 3,4,5-Trichlorophenol already formulated into the standard)

Vial 1 (Methanol)**Labeled Compounds****Concentration (µg/mL)**2,4-Dichlorophenol (¹³C₆,99%) 2504-Chloroguaiacol (¹³C₆,99%) 2504,5-Dichlorocatechol (Ring-¹³C₆,99%) 2504,5,6-Trichloroguaiacol (¹³C₆,99%) 250Pentachlorophenol (¹³C₆,99%) 2503,4,5,6-Tetrachloroguaiacol (¹³C₆,99%) 2503,4,5,6-Tetrachlorocatechol (Ring-¹³C₆,99%) 250

3,4,5-Trichlorophenol (Unlabeled) 250

Vial 2 (Acetone)**Labeled Compound****Concentration (µg/mL)**5-Chlorovanillin (Ring-¹³C₆,99%) 250