

ES-5438

PAH Surrogates Standard Mixture Unlabeled

10 mL
90% Toluene/
10% Isooctane

Unlabeled Component	Concentration (ug/mL)
Naphthalene	200
Benz[<i>a</i>]anthracene	200
Phenanthrene	200
Fluoranthene	200
Benzo[<i>b</i>]fluoranthene	200
Benzo[<i>a</i>]pyrene	200
Benzo[<i>g,h,i</i>]perylene	200
Indeno[1,2,3- <i>cd</i>]pyrene	200
Dibenz[<i>a,h</i>]anthracene	200
Acenaphthylene	200
Acenaphthene	200
Fluorene	200
Pyrene	200
Benzo[<i>k</i>]fluoranthene	200
Perylene	200
Chrysene	200