



Cambridge Isotope Laboratories, Inc.

3 Highwood Drive, Tewksbury, MA 01876 USA

PH: 800.322.1174 (N. America) PH: 978-749-8000

Fax: 978-749-2768 WEB: www.isotope.com

ES-5630

JECS Phenol/Paraben Syringe Standard

1.2 mL in methanol-OD

| Labeled | Concentration (ng/mL) |
|---|----------------------------------|
| Methyl paraben (methyl 4-hydroxybenzoate) (2,3,5,6-D ₄ , 98%) | 250 |
| <i>n</i> -Propyl paraben (<i>n</i> -propyl 4-hydroxybenzoate) (2,3,5,6-D ₄ , 98%) | 250 |
| <i>p</i> - <i>n</i> -Nonylphenol (ring-D ₄ , OD, 98%) | 2500 |
| Bisphenol S (ring-D ₈ , 98%) | 250 |
| Bisphenol F (D ₁₀ , 98%) | 2500 |
| Bisphenol A (3,3',5,5'-D ₄ , 97%) | 2500 |